

Anti-NGFR/p75 neurotrophin receptor antibody [MC192] (STJA0040799)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	In-vivo
Host / Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified
Formulation	Lyophilized from a 1 mg/mL solution containing PBS pH 7.2-7.6 without preservative.
Isotype	IgG1
Storage Instruction	Spin vial briefly before opening. Reconstitute in 50 μ L sterile-filtered, ultrapure water to achieve a conjugate concentration of 1 mg/mL. Centrifuge to remove any insoluble material. Final buffer contains no preservatives. Lyophilized product is ship

TARGET INFORMATION

Gene ID	24596
Gene Symbol	Ngfr
UniProt ID	TNR16_RAT
Immunogen Sequence	Rat
Specificity	MC192 antibody is specific only for rat NGFR, no reactivity to human or mouse NGFR has been reported This monoclonal antibody does not cross react with p75NTR-expressing cells in other species than rat.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---