

Anti-S-arrestin antibody (C-Term) [S128] (STJA0040804)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Mouse
Reactivity	Bovine/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 3.5138888888888888 IHC 0.73611111111111116
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide.
Isotype	IgG1
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution divide into aliquot

TARGET INFORMATION

Immunogen Region	C-Term
Immunogen Sequence	Bovine
Specificity	Bovine, reacts with pig. Other species not yet tested. Antibody detects arrestin-1 protein only and does not cross-react with the other 3 arrestin molecules.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---