

Anti-Methyl-CpG binding protein 2 antibody (Full-Length) [4F11] (STJA0040809)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:5000-1:10000 IHC 1:1000-1:5000
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide.
Isotype	IgG1
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution divide into aliquo

TARGET INFORMATION

Gene ID	4204
Gene Symbol	MECP2
UniProt ID	MECP2_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Reacts with human, mouse, rat.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---