

Anti-Rhodopsin antibody [B630] (STJA0040811)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Mouse
Reactivity	Bovine/Human/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 3.5138888888888888 IHC 0.73611111111111116
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide.
Isotype	IgG1
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution divide into aliquo

TARGET INFORMATION

Immunogen Sequence	Bovine
Specificity	Bovine, reacts with human, rat, mouse, cow, pig. Expected to react with other mammalian species.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---