

Anti-Beta-synuclein antibody (C-Term) [6A10] (STJA0040819)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Mouse
Reactivity	Bovine/Human/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 0.73611111111111116 IHC 1:500-1:1000
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG1
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water to achieve a 1 mg/mL concentration. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution divide into aliquo

TARGET INFORMATION

Gene ID	6620
Gene Symbol	SNCB
UniProt ID	SYUB_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	Validated by Western blotting and immunohistochemical procedures. Specific for Beta-synuclein, does not cross-react with Alpha-or Gamma-synuclein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---