

Anti-Choline O-acetyltransferase antibody (aa: 167-188) (STJA0040829)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/ICC/IHC-F/IHC-P
Host / Source	Rabbit
Reactivity	Guinea Pig/Mouse/Pig/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 1:400-1:2000 ICC 1:400-1:2000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Immunogen Region	aa: 167-188
Immunogen Sequence	Pig
Specificity	This antiserum stains cholinergic neurons in guinea-pig and rabbit. This antiserum is known to react with ChAT of origin guinea-pig, mouse, rat and rabbit.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------