

## Anti-Capsaicin receptor antibody (aa: 531-541) (STJA0040830)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 1:1000-1:2000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 $\mu$ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

### TARGET INFORMATION

Gene ID	<a href="#">7442</a>
Gene Symbol	<a href="#">TRPV1</a>
UniProt ID	<a href="#">TRPV1_HUMAN</a>
Immunogen Region	aa: 531-541
Immunogen Sequence	Human
Specificity	Immunohistochemical analysis of rat DRG indicates a superb degree of specificity for this serum. This antibody is known to react with rat Capsaicin and due to sequence homology is expected to react with mouse Capsaicin. Other species have not yet been

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---