

## Anti-Gamma-synuclein antibody (aa: 114-127) (STJA0040832)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Rabbit
Reactivity	Human/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:500-1:1000 IHC 1:500-1:1000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 $\mu$ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

### TARGET INFORMATION

Gene ID	<a href="#">6623</a>
Gene Symbol	<a href="#">SNCG</a>
UniProt ID	<a href="#">SYUG_HUMAN</a>
Immunogen Region	aa: 114-127
Immunogen Sequence	Human
Specificity	Immunohistochemical/Western blot analysis indicate a high level of specificity for this antiserum for gamma synuclein. This antiserum is known to react with human and rat gamma synuclein.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---