

Anti-Glial cell-derived neurotrophic factor receptor alpha-2 antibody (aa: 376-390) (STJA0040835)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:500-1:1000 IHC 1:500-1:1000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	2675
Gene Symbol	GFRA2
UniProt ID	GFRA2_HUMAN
Immunogen Region	aa: 376-390
Immunogen Sequence	Human
Specificity	Western blot analysis indicates a high level of specificity to GFR alpha-2 for this antiserum. This antibody is known to react with rat GFR alpha-2 and from IHC analysis, there appears to be no cross-reactivity with other members of the GFR family.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---