

Anti-Adenosine triphosphate-ase antibody (aa: 147-161) (STJA0040840)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 1:500-1:2000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	74772
Gene Symbol	Atp13a2
UniProt ID	AT132_MOUSE
Immunogen Region	aa: 147-161
Immunogen Sequence	Mouse
Specificity	Immunohistochemistry shows specific staining for ATPase. This antiserum is known to react with rat ATPase. Other species have not yet tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---