

Anti-Contactin-6 antibody (GYSPHYRWKQ) (STJA0040863)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 1:1000-1:4000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for increased stability, and at 2-8°C with an appropriate antibacteria

TARGET INFORMATION

Gene ID	53870
Gene Symbol	Cntn6
UniProt ID	CNTN6_MOUSE
Immunogen Region	GYSPHYRWKQ
Immunogen Sequence	Mouse
Specificity	This antiserum stains beautifully the NB3 in DRG, where it nicely stains a subpopulation of neurons and in the trigeminal nerve root, where staining is restricted to the oligodendrocytes. Confirmed to react with rat Contactin-6 and expected to recogn

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---