

Anti-Parkinson disease protein 2 antibody (RILGEEQYNRYQQYGAEEC) (STJA0040868)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IHC-F/WB
Host / Source	Rabbit
Reactivity	Guinea Pig/Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:500-1:2000 IHC 1:500-1:2000 ICC 1:500-1:2000
Isotype	Mixed
Molecular Weight	50 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-

TARGET INFORMATION

Gene ID	5071
Gene Symbol	PRKN
UniProt ID	PRKN_HUMAN
Immunogen Region	RILGEEQYNRYQQYGAEEC
Immunogen Sequence	Human
Specificity	This antiserum is known to specifically recognize Parkin shown by IHC and WB. This antibody is known to react with Parkin of guinea pig and rat.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---