

**Anti-Noxa antibody (aa: 1-16) (STJA0040876)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Rabbit             |
| Reactivity    | Mouse              |

**PRODUCT PROPERTIES**

|                     |                                                                                                                                                                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                                                                                                                       |
| Conjugation         | Unconjugated                                                                                                                                                                                                                                                     |
| Purification        | Whole serum                                                                                                                                                                                                                                                      |
| Dilution Range      | WB 1:1000-1:2000                                                                                                                                                                                                                                                 |
| Isotype             | Mixed                                                                                                                                                                                                                                                            |
| Storage Instruction | Spin vial briefly before opening. Reconstitute in 100 $\mu$ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial |

**TARGET INFORMATION**

|                    |                                                                                                                                                                                                           |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">58801</a>                                                                                                                                                                                     |
| Gene Symbol        | <a href="#">Pmaip1</a>                                                                                                                                                                                    |
| UniProt ID         | <a href="#">APR_MOUSE</a>                                                                                                                                                                                 |
| Immunogen Region   | aa: 1-16                                                                                                                                                                                                  |
| Immunogen Sequence | Mouse                                                                                                                                                                                                     |
| Specificity        | Western blot analysis of cells infected with Noxa adenoviruses and BAF indicates a high level of specificity for this antiserum. This antiserum cross-reacts with mouse. Not yet tested in other species. |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|