

Anti-Noxa antibody (aa: 75-90) (STJA0040877)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:1000-1:2000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	58801
Gene Symbol	Pmaip1
UniProt ID	APR_MOUSE
Immunogen Region	aa: 75-90
Immunogen Sequence	Mouse
Specificity	Western blot analysis of cells infected with Noxa adenoviruses and BAF indicates a high level of specificity for this antiserum. This antiserum cross-reacts with mouse. Not yet tested in other species.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---