

Anti-Superoxide dismutase 1 antibody (aa: 24-36) (STJA0040879)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:500-1:1000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	20655
Gene Symbol	Sod1
UniProt ID	SODC_MOUSE
Immunogen Region	aa: 24-36
Immunogen Sequence	Mouse
Specificity	Specificity was confirmed by Western blot detecting mouse superoxide dismutase (SOD1). This antiserum is known to react with mouse and rat superoxide dismutase (SOD1) protein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---