

**Anti-Contactin-6 antibody (KLHKMSSLSSVGVQILKPSTQF) (STJA0040880)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Rabbit
Reactivity	Rat

**PRODUCT PROPERTIES**

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 1:1000-1:5000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 $\mu$ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at minus 20°C for increased stability, and at 2-8°C with an appropriate antiba

**TARGET INFORMATION**

Gene ID	<a href="#">27256</a>
Gene Symbol	<a href="#">Cntn6</a>
UniProt ID	<a href="#">CNTN6_RAT</a>
Immunogen Region	KLHKMSSLSSVGVQILKPSTQF
Immunogen Sequence	Rat
Specificity	This antiserum stains beautifully the NB3 in DRG, where it nicely stains a subpopulation of neurons and in the trigeminal nerve root, where staining is restricted to the oligodendrocytes. Confirmed to react with rat Contactin-6.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---