

## Anti-D Tag antibody (DYKDDDDK) (STJA0040881)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Rabbit
Reactivity	Species Independent

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 0.7361111111111116 IHC 0.1805555555555555 ICC 0.1805555555555555 IP 0.1805555555555555
Formulation	PBS, pH 7.4 with 0.05% sodium azide
Isotype	IgG
Storage Instruction	Stable for 1 year at -20°C from the date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.

### TARGET INFORMATION

Immunogen Region	DYKDDDDK
Immunogen Sequence	Species Independent
Specificity	Anti D-tag antibody recognizes over-expressed proteins containing the D-tag fused to either amino- or carboxy-termini of targeted proteins in transfected mammalian cells.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---