

Anti-Adiponectin Receptor 1 antibody (4-142 aa) (STJA0040886)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1.4305555555555556
Formulation	15% glycerol
Isotype	Mixed
Storage Instruction	Aliquot and store at -20°C for a higher stability. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	72674
Gene Symbol	Adipor1
UniProt ID	PAQR1_MOUSE
Immunogen Region	4-142 aa
Immunogen Sequence	Mouse
Specificity	This antiserum is known to recognize both mouse and human Adiponectin Receptor 1.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---