

Anti-Adiponectin Receptor 2 antibody (78-219 aa) (STJA0040887)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1.4305555555555556
Formulation	15% glycerol
Isotype	Mixed
Storage Instruction	Aliquot and store at -20°C for a higher stability. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79602
Gene Symbol	ADIPOR2
UniProt ID	PAQR2_HUMAN
Immunogen Region	78-219 aa
Immunogen Sequence	Human
Specificity	This antiserum is known to recognize both mouse and human Adiponectin Receptor 2.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---