

Anti-Adaptor protein, phosphotyrosine interacting with PH domain and leucine zipper 1 antibody (C-Term) (STJA0040888)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 0.7361111111111116 IP 3 µg/reaction
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. At least 12 months after purchase at 2-8°C (lyophilized formulations). After reconstitution, aliquot and store at

TARGET INFORMATION

Gene ID	26060
Gene Symbol	APPL1
UniProt ID	DP13A_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	This antiserum is known to recognize both mouse and human APPL1.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---