

## Anti-Coronin-1A antibody (C-Term) (STJA0040891)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Rabbit
Reactivity	Bovine/Human/Mouse/Pig/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:2500-1:5000 IP 1:500-1:1000
Formulation	Lyophilized with sodium azide.
Isotype	Mixed
Molecular Weight	51 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 50 $\mu$ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for a higher stability. Avoid freeze-thaw.

### TARGET INFORMATION

Gene ID	<a href="#">11151</a>
Gene Symbol	<a href="#">CORO1A</a>
UniProt ID	<a href="#">COR1A_HUMAN</a>
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	The specificity of this antibody has been confirmed by WB and ICC against the antigen. Human, rat, mouse and feline. Predicted to react with other mammalian tissue.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---