

Anti-Autophagy-related protein 12 antibody (N-Term) (STJA0040892)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 1-5 µg/mL
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute with 500 µL of sterile-filtered, ultra pure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibac

TARGET INFORMATION

Gene ID	9140
Gene Symbol	ATG12
UniProt ID	ATG12_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	IHC and WB confirmed the specificity for ATG12. Human, rat.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---