

Anti-Internexin alpha antibody (STJA0040894)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Rabbit
Reactivity	Cat/Human/Mouse/Other Mammals/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:10000-1:20000 IHC 1:500-1:1000 ICC 1:500-1:1000
Formulation	Lyophilized with sodium azide.
Isotype	Mixed
Molecular Weight	64-66 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 50 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for a higher stability. Avoid freeze-tha

TARGET INFORMATION

Gene Symbol	Ina Uniprot ID= AINX_RAT "
Immunogen Sequence	Rat
Specificity	The specificity of this antibody has been confirmed by WB. This antibody is specific for the 64-66 kDa alpha-internexin protein. Molecular weight will depend on species. Hu, rat, Ms, Fel, and other mammals.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---