

Anti-Neurofilament light polypeptide antibody (Full-Length) (STJA0040899)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ICC/IHC/WB
Host / Source	Rabbit
Reactivity	Bovine/Human/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:10000-1:15000 IHC 3.5138888888888888 ICC 1:2000-1:5000
Formulation	Lyophilized with sodium azide.
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute with 50 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C After reconstitution of lyophilized antibody, aliquot and store at -20°C for

TARGET INFORMATION

Gene ID	4747
Gene Symbol	NEFL
UniProt ID	NFL_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Species cross-reactivity includes human, rat, mouse, cow, and pig.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---