

Anti-Microtubule-associated protein 1 light chain 3 gamma antibody (N-Term) (STJA0040902)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:100-1:1000 IHC 1:100-1:1000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	440738
Gene Symbol	MAP1LC3C
UniProt ID	MLP3C_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	IHC and WB confirmed the specificity for MAP1LC3 C. This antibody should recognize MAP1LC3 C only and not the other forms MAP1LC3A and MAP1LC3B. Human, other species have not yet been tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---