

Anti-Peripherin antibody (Full-Length) (STJA0040903)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Rabbit
Reactivity	Bovine/Human/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:2000-1:10000 ICC 1:1000-1:2000
Formulation	Lyophilized with sodium azide.
Isotype	Mixed
Molecular Weight	48 kDa, 57 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute with 50 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for a higher stability. Avoid freeze-thaw.

TARGET INFORMATION

Gene Symbol	Prph Uniprot ID= PERL_RAT "
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	The specificity of this antibody has been confirmed by WB. This antibody detects ~57 kDa Peripherin protein and a smaller molecule derived from Peripherin at ~48 kDa. Hu, rat, Ms, Fel, Rb. Predicted to react with other mammalian tissue.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---