

Anti-Gamma-aminobutyric acid receptor-associated protein-like 2 antibody (N-Term) (STJA0040907)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:200-1:1000 IHC 1:200-1:1000 ICC 1:200-1:1000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	11345
Gene Symbol	GABARAPL2
UniProt ID	GBRL2_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	IHC, WB and ELISA confirmed the specificity for GABARAP L2 (GABARAPL2). A 14 kDa band, that corresponds to the molecular weight of GABARAPL2, is detected via Western blot analysis. Human, rat. Other species not yet tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---