

Anti-Gamma-aminobutyric acid receptor-associated protein antibody (N-Term) (STJA0040908)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:200-1:1000 IHC 1:200-1:1000 ICC 1:200-1:1000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	11337
Gene Symbol	GABARAP
UniProt ID	GBRAP_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	IHC, WB and ELISA confirmed the specificity for GABARAP. Human, rat. Other species not yet tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---