

Anti-Autophagy-related protein 10 antibody (C-Term) (STJA0040909)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 1:100-1:1000 ICC 1:100-1:1000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at minus 20°C for a higher stability, and at 2-8°C with an appropriate antibac

TARGET INFORMATION

Gene ID	83734
Gene Symbol	ATG10
UniProt ID	ATG10_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	IHC, WB and ELISA confirmed the specificity for ATG10. Human. Other species not yet tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---