

Anti-Microtubule-associated proteins 1A/1B light chain 3A antibody (N-Term) (STJA0040910)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 0.388888888888889 IHC 1:100-1:1000 ICC 1:100-1:1000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at-20°C for higher stability, and at 2-8°C with an appropriate antibacterial a

TARGET INFORMATION

Gene ID	84557
Gene Symbol	MAP1LC3A
UniProt ID	MLP3A_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	IHC, WB and ELISA confirmed the specificity for MAP1LC3 A. Human, rat. Other species not yet tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---