

Anti-Adrenocorticotrop hormone antibody (STJA0040913)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Rabbit
Reactivity	Human/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	IHC 5-15 μ g/mL ICC 5-15 μ g/mL
Formulation	Lyophilized with BSA.
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 0.05 mL of sterile-filtered PBS (pH 7.4) to achieve an antibody concentration of 1 mg/mL. Centrifuge to remove any insoluble material. At least 12 months after purchase at 2-8°C (lyophilized formulati

TARGET INFORMATION

Immunogen Sequence	Pig
Specificity	Human; mouse; rat; ACTH is highly conserved so cross-reactivity with other species is expected.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---