

Anti-Sortilin antibody (STJA0040922)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC-F
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 0.73611111111111116
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at 20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	6272
Gene Symbol	SORT1
UniProt ID	SORT_HUMAN
Immunogen Sequence	Human
Specificity	Specificity was demonstrated by immunohistochemistry. Reacts with human, rat and mouse Sortilin. Other species not yet tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---