

Anti-Somatostatin antibody (103-116) (STJA0040924)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/IHC-P
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	IHC 5-15 μ g/mL ICC 5-15 μ g/mL
Formulation	Lyophilized with BSA.
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 0.05 mL of sterile-filtered PBS (pH 7.4). Centrifuge to remove any insoluble material. At least 12 months after purchase at 2-8°C (lyophilized formulations). After reconstitution, aliquot and store at

TARGET INFORMATION

Gene ID	6750
Gene Symbol	SST
UniProt ID	SMS_HUMAN
Immunogen Region	103-116
Immunogen Sequence	Human
Specificity	Human Somatostatin. Cross reactivity is expected with mouse, rat, other species not yet tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---