

Anti-Tyrosine Kinase Receptor C antibody (STJA0040926)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 1:1000-1:3000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	4916
Gene Symbol	NTRK3
UniProt ID	NTRK3_HUMAN
Immunogen Sequence	Human
Specificity	Specificity was demonstrated by immunohistochemistry. This antibody was used to stain cryostat sections of the rat peripheral sensory ganglia. Reacts with human, rat and mouse TrkC. Other species have not yet been tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---