

Anti-Histone H3 antibody (C-Term) (STJA0040930)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:1000-1:5000
Formulation	PBS with 0.02% sodium azide.
Isotype	IgG
Storage Instruction	Short term storage at 2-8°C for one week. At-20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	8350/8351/8352/8353/8354/8355/8356/8357/8358/8968
Gene Symbol	H3C1.H3C2.H3C3.H3C4.H3C6.H3C7.H3C8.H3C10.H3C11.H3C12
UniProt ID	H31_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	Antibody binds to various Histone H3 proteins. Human, mouse, rat, bovine.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---