

## Anti-Leptin antibody (STJA0040933)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:1000-1:2000
Formulation	PBS with 0.02% sodium azide.
Isotype	Mixed
Storage Instruction	Short term storage at 2-8°C for one week. At-20°C as an undiluted liquid for up to 12 months.

### TARGET INFORMATION

Gene ID	<a href="#">3952</a>
Gene Symbol	<a href="#">LEP</a>
UniProt ID	<a href="#">LEP_HUMAN</a>
Immunogen Sequence	Human
Specificity	Human Leptin.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---