

Anti-Neurturin antibody (STJA0040937)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1.4305555555555556
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	4902
Gene Symbol	NRTN
UniProt ID	NRTN_HUMAN
Immunogen Sequence	Human
Specificity	Specificity of this product was demonstrated by immunohistochemistry and dot blot showed no cross reactivity with GDNF. This antibody is known to react with human, mouse and rat.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---