

Anti-Autophagy-related protein 16-1 antibody (N-Term) (STJA0040940)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/ICC/WB
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:1000-1:4000 IHC 1:200-1:1000
Isotype	Mixed
Molecular Weight	22 kDa
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	55054
Gene Symbol	ATG16L1
UniProt ID	A16L1_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	Detects human protein via non-reduced WB, and ELISA human, rat. Other species not yet tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------