

Anti-Glutathione peroxidase 1 antibody (GPx-1) (STJA0040941)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC-P/WB
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:1000-1:2000 IHC 1:200-1:1000
Formulation	Lyophilized serum with 0.02% thimerosal
Isotype	Mixed
Molecular Weight	22 kDa (see application details)
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. At least 12 months after purchase at 2-8°C (lyophilized formulations). After reconstitution, aliquot and store at

TARGET INFORMATION

Gene ID	2876
Gene Symbol	GPX1
UniProt ID	GPX1_HUMAN
Immunogen Region	GPx-1
Immunogen Sequence	Human
Specificity	Human and rat cellular selenium dependent glutathione peroxidase. Human; rat,;

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---