

Anti-Sequestosome 1 antibody (C-Term) (STJA0040949)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Rabbit
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	IHC 1:500-1:5000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at -8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	8878
Gene Symbol	SQSTM1
UniProt ID	SQSTM_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	The specificity was confirmed by immunohistochemical analysis on paraffin embedded human brain. This antibody is expected to recognize p60 in mouse, rat and bovine.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---