

Anti-AP-1 complex subunit beta-1 antibody (11-29 aa) (STJA0040955)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ICC/IHC-F/IHC-P/WB
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:100-1:10000 IHC 1:1000-1:5000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	162
Gene Symbol	AP1B1
UniProt ID	AP1B1_HUMAN
Immunogen Region	11-29 aa
Immunogen Sequence	Human
Specificity	This antibody has been shown to be specific for AP-1 complex subunit beta-1. Rat and human, other species have not yet tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---