

Anti-Conserved dopamine neurotrophic factor antibody (STJA0040956)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:6000-1:10000
Formulation	Concentrated ammonium sulphate in PBS pH 7.4 containing no preservatives.
Isotype	IgG
Storage Instruction	Store antibody stock solution at +4°C upon receipt. DO NOT FREEZE. This product is an (NH ₄) ₂ SO ₄ precipitate. To use, mix well by pipetting or vortexing gently prior use, then withdraw the desired amount of antibody stock solution and dilute into you

TARGET INFORMATION

Gene ID	441549
Gene Symbol	CDNF
UniProt ID	CDNF_HUMAN
Immunogen Sequence	Human
Specificity	CDNF (conserved dopamine neurotrophic factor).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---