

Anti-Superoxide dismutase 1 antibody (SOD1) (STJA0040958)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/IHC-P/WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:1000-1:4000 IHC 1:100-1:1000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	6647
Gene Symbol	SOD1
UniProt ID	SODC_HUMAN
Immunogen Region	SOD1
Immunogen Sequence	Human
Specificity	This antibody has been shown to be specific for superoxide dismutase (SOD1) protein. Human. Does not detect mouse and rat SOD1.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---