

Anti-Presenilin 1 loop region antibody (C-Term) (STJA0040959)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 0.7361111111111116
Formulation	Lyophilized from PBS, pH 7.4. Contains no preservative.
Isotype	IgG
Molecular Weight	19 kDa, 45 kDa (See application details.)
Storage Instruction	Spin vial briefly before opening. Reconstitute in 500 μ L sterile-filtered 1X PBS, pH 7.2-7.6. Centrifuge to remove any insoluble material. Short term storage at 2-8°C for one week. At -20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	5663
Gene Symbol	PSEN1
UniProt ID	PSN1_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	Confirmed by Western blotting using transfected cells, presenilin 1 knock-out mouse cells and mouse and human brain.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---