

Anti-Nicastrin, N-terminal antibody (N-Term) (STJA0040965)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	WB 3-10 µg/mL
Formulation	Lyophilized from PBS, pH 7.4. Contains no preservative.
Isotype	IgG
Molecular Weight	145 kDa (See application details.)
Storage Instruction	Spin vial briefly before opening. Reconstitute in 500 µL sterile-filtered 1X PBS, pH 7.2-7.6. Centrifuge to remove any insoluble material. Short term storage at 2-8°C for one week. At 20°C as an undiluted liquid for up to 12 months.

TARGET INFORMATION

Gene ID	23385
Gene Symbol	NCSTN
UniProt ID	NICA_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	Confirmed by WB using peptide absorption.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---