

## Anti-Peroxiredoxin-6 antibody (STJA0040967)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/ICC/IHC-F/IHC-P/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:100-1:1000 IHC 1:1000-1:4000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 $\mu$ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

### TARGET INFORMATION

Gene ID	<a href="#">94167</a>
Gene Symbol	<a href="#">Prdx6</a>
UniProt ID	<a href="#">PRDX6_RAT</a>
Immunogen Sequence	Human
Specificity	This antibody has been shown to be specific for Peroxiredoxin-6 protein. Rat, human and mouse, other species have not yet been tested.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---