

Anti-Neuropeptide Y antibody (29-64 aa) (STJA0040975)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Sheep

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	IHC 1:2000-1:4000
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 50 μ L sterile-filtered, ultrapure water. Pulse-centrifuge to remove any insoluble material. Maintain lyophilized material at 2-8°C for up to 12 months. After reconstitution keep aliquots at -20°C for u

TARGET INFORMATION

Immunogen Region	29-64 aa
Immunogen Sequence	Sheep
Specificity	Specificity has been confirmed by direct ELISA against the antigen. Sheep, human, mouse & rat.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---