

## Anti-Activity-regulated gene 3.1 protein homolog antibody (STJA0040978)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	IHC 2-4 µg/mL
Formulation	Lyophilized from a solution containing PBS pH 7.4
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 500 µL sterile-filtered 1X PBS, pH 7.2-7.6. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacteri

### TARGET INFORMATION

Gene ID	<a href="#">23237</a>
Gene Symbol	<a href="#">ARC</a>
UniProt ID	<a href="#">ARC_HUMAN</a>
Immunogen Sequence	Human
Specificity	Specific for Arg3.1. Rat and mouse. Other species not yet tested.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---