

Anti-Beta casein protein antibody (STJA0040982)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/WB |
| Host / Source | Rabbit |
| Reactivity | Bovine |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Polyclonal |
| Conjugation | Unconjugated |
| Purification | Affinity-purified |
| Formulation | PBS, pH 7.4 |
| Isotype | IgG |
| Storage Instruction | Spin vial briefly before opening and centrifuge to remove any insoluble material. After opening maintain at -20°C in undiluted aliquots for up to 6 months. For short-term storage, keep aliquot at 2-8°C for up to one week. Avoid repeated freeze-thaw cy |

TARGET INFORMATION

| | |
|--------------------|--|
| Immunogen Sequence | Bovine |
| Specificity | Bovine beta-casein. Detects both A1 and A2 beta-casein isoforms. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|