

Anti-2 NA, 3 NA-cyclic-nucleotide 3 NA-phosphodiesterase antibody (Full-Length) (STJA0040984)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:1000-1:5000 ICC 1:1000-1:2000
Formulation	Lyophilized with sodium azide.
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution divide into aliquots and store at 20°C for long-term

TARGET INFORMATION

Gene ID	1267
Gene Symbol	CNP
UniProt ID	CN37_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Reacts with human, rat, mouse.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---